

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO
Test: Acid Phosphatase
Catalog # : A7503

UNITS	U/L
CALCULATION FACTOR	969
STANDARD 1 CONC.	0
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	0.4
TEMPERATURE (DEG C)	30
TYPE OF ANALYSIS	2
WAVELENGTH (NM)	405
SAMPLE VOLUME (UL)	20
DILUENT VOLUME (UL)	50
REAGENT VOLUME (UL)	150
INCUBATION TIME (SEC)	0
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	300
TIME INTERVAL (SEC)	20
NUMBER OF READINGS	15
BLANKING MODE	0
PRINTOUT MODE	1
ALPHA CODE	11

Reconstitute Acid phosphatase vials 7.5 ml of distilled water. Swirl to dissolve.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO

Test: Albumin

Catalog # : A7502

UNITS	G/DL
CALCULATION FACTOR	0
STANDARD 1 CONC.	*
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	6
TEMPERATURE (DEG C)	25
TYPE OF ANALYSIS	1
WAVELENGTH (NM)	630
SAMPLE VOLUME (UL)	3
DILUENT VOLUME (UL)	15
REAGENT VOLUME (UL)	300
INCUBATION TIME (SEC)	10
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	0.5
TIME INTERVAL (SEC)	10
NUMBER OF READINGS	1
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

Reagent supplied as a ready to use liquid

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO
Test: Alkaline Phosphatase
Catalog # : A7505/A7516

UNITS	U/L
CALCULATION FACTOR	2660
STANDARD 1 CONC.	0
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	0.5
TEMPERATURE (DEG C)	37
TYPE OF ANALYSIS	2
WAVELENGTH (NM)	405
SAMPLE VOLUME (UL)	5
DILUENT VOLUME (UL)	30
REAGENT VOLUME (UL)	150
INCUBATION TIME (SEC)	0
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	50
TIME INTERVAL (SEC)	10
NUMBER OF READINGS	16
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

Add 13 ml distilled water to the 15 ml size vial and 43 ml to the 50 ml size vial.
It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO

Test: ALT

Catalog # : A7525/A7526

UNITS	U/L
CALCULATION FACTOR	2010
STANDARD 1 CONC.	0
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	0.3
TEMPERATURE (DEG C)	37
TYPE OF ANALYSIS	2
WAVELENGTH (NM)	340
SAMPLE VOLUME (UL)	20
DILUENT VOLUME (UL)	50
REAGENT VOLUME (UL)	150
INCUBATION TIME (SEC)	0
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	50
TIME INTERVAL (SEC)	10
NUMBER OF READINGS	16
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	13

Add 13 ml distilled water to the 15 ml size vial and 43 ml to the 50 ml size vial.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO

Test: Amylase

Catalog # : A7564

UNITS	U/L
CALCULATION FACTOR	2769
STANDARD 1 CONC.	0
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	2
TEMPERATURE (DEG C)	37
TYPE OF ANALYSIS	2
WAVELENGTH (NM)	405
SAMPLE VOLUME (UL)	7
DILUENT VOLUME (UL)	15
REAGENT VOLUME (UL)	265
INCUBATION TIME (SEC)	0
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	15
TIME INTERVAL (SEC)	10
NUMBER OF READINGS	10
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

Reagent supplied as a ready to use liquid

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO

Test: AST

Catalog # : A7560/A7561

UNITS	U/L
CALCULATION FACTOR	2010
STANDARD 1 CONC.	0
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	0.3
TEMPERATURE (DEG C)	37
TYPE OF ANALYSIS	2
WAVELENGTH (NM)	340
SAMPLE VOLUME (UL)	20
DILUENT VOLUME (UL)	50
REAGENT VOLUME (UL)	150
INCUBATION TIME (SEC)	0
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	50
TIME INTERVAL (SEC)	10
NUMBER OF READINGS	16
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	13

Add 13 ml distilled water to the 15 ml size vial and 43 ml to the 50 ml size vial.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO

Test: Direct Bilirubin

Catalog # : B7538

UNITS	MG/DL
CALCULATION FACTOR	0
STANDARD 1 CONC.	*
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	18
TEMPERATURE (DEG C)	25
TYPE OF ANALYSIS	6
WAVELENGTH (NM)	560
SAMPLE VOLUME (UL)	20
DILUENT VOLUME (UL)	50
REAGENT VOLUME (UL)	200
INCUBATION TIME (SEC)	10
START REAGENT VOLUME	2
TIME OF FIRST READING (SEC)	0.5
TIME INTERVAL (SEC)	300
NUMBER OF READINGS	2
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

Place Direct Bilirubin in reagent tray with nitrite reagent in the second reagent cup

Do not prepare a working reagent.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO
Test: Total Bilirubin
Catalog # : B7576

UNITS	MG/DL
CALCULATION FACTOR	0
STANDARD 1 CONC.	*
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	18
TEMPERATURE (DEG C)	37
TYPE OF ANALYSIS	6
WAVELENGTH (NM)	560
SAMPLE VOLUME (UL)	15
DILUENT VOLUME (UL)	10
REAGENT VOLUME (UL)	290
INCUBATION TIME (SEC)	10
START REAGENT VOLUME	3
TIME OF FIRST READING (SEC)	0.5
TIME INTERVAL (SEC)	60
NUMBER OF READINGS	2
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

Place Total Bilirubin in reagent tray with nitrite reagent in the second reagent cup
Do not prepare a working reagent.
It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO
Test: BUN
Catalog # : B7552

UNITS	MG/DL
CALCULATION FACTOR	0
STANDARD 1 CONC.	*
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	80
TEMPERATURE (DEG C)	30
TYPE OF ANALYSIS	4
WAVELENGTH (NM)	340
SAMPLE VOLUME (UL)	3
DILUENT VOLUME (UL)	30
REAGENT VOLUME (UL)	300
INCUBATION TIME (SEC)	0
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	30
TIME INTERVAL (SEC)	10
NUMBER OF READINGS	4
BLANKING MODE	1
PRINTOUT MODE	2
ALPHA CODE	13

Add 13.5 ml distilled water to the 15 ml size vial and 45 ml to the 50 ml size vial.
It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO
Test: Calcium (Arsenazo)
Catalog # : C7529

UNITS	MG/DL
CALCULATION FACTOR	0
STANDARD 1 CONC.	10
STANDARD 2 CONC.	10
STANDARD 3 CONC.	10
LIMIT	15
TEMPERATURE (DEG C)	37
TYPE OF ANALYSIS	1
WAVELENGTH (NM)	630
SAMPLE VOLUME (UL)	2
DILUENT VOLUME (UL)	15
REAGENT VOLUME (UL)	200
INCUBATION TIME (SEC)	10
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	1
TIME INTERVAL (SEC)	60
NUMBER OF READINGS	2
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

Reagent supplied as a ready to use liquid

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO
Test: Chloride
Catalog # : C7501

UNITS	mEq/L
CALCULATION FACTOR	0
STANDARD 1 CONC.	*
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	130
TEMPERATURE (DEG C)	30
TYPE OF ANALYSIS	1
WAVELENGTH (NM)	480
SAMPLE VOLUME (UL)	2
DILUENT VOLUME (UL)	10
REAGENT VOLUME (UL)	300
INCUBATION TIME (SEC)	10
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	0.5
TIME INTERVAL (SEC)	30
NUMBER OF READINGS	2
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

Prepare reagent according to package insert.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO

Test: Cholesterol

Catalog # : C7509

UNITS	MG/DL
CALCULATION FACTOR	0
STANDARD 1 CONC.	*
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	600
TEMPERATURE (DEG C)	37
TYPE OF ANALYSIS	5
WAVELENGTH (NM)	520
SAMPLE VOLUME (UL)	4
DILUENT VOLUME (UL)	30
REAGENT VOLUME (UL)	250
INCUBATION TIME (SEC)	0
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	0.5
TIME INTERVAL (SEC)	600
NUMBER OF READINGS	2
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO
Test: CPK
Catalog # : C7512

UNITS	U/L
CALCULATION FACTOR	8039
STANDARD 1 CONC.	0
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	0.3
TEMPERATURE (DEG C)	37
TYPE OF ANALYSIS	2
WAVELENGTH (NM)	340
SAMPLE VOLUME (UL)	5
DILUENT VOLUME (UL)	30
REAGENT VOLUME (UL)	200
INCUBATION TIME (SEC)	0
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	180
TIME INTERVAL (SEC)	10
NUMBER OF READINGS	10
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

Add 13 ml distilled water to the 15 ml size vial and 43 ml to the 50 ml size vial.
It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO
Test: Creatinine
Catalog # : C7539

UNITS	MG/DL
CALCULATION FACTOR	0
STANDARD 1 CONC.	*
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	25
TEMPERATURE (DEG C)	30
TYPE OF ANALYSIS	5
WAVELENGTH (NM)	510
SAMPLE VOLUME (UL)	10
DILUENT VOLUME (UL)	30
REAGENT VOLUME (UL)	250
INCUBATION TIME (SEC)	0
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	20
TIME INTERVAL (SEC)	60
NUMBER OF READINGS	2
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

Prepare reagent according to package insert.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO
Test: GGT
Catalog # : G7570/G7571

UNITS	U/L
CALCULATION FACTOR	2631
STANDARD 1 CONC.	0
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	0.4
TEMPERATURE (DEG C)	37
TYPE OF ANALYSIS	2
WAVELENGTH (NM)	405
SAMPLE VOLUME (UL)	10
DILUENT VOLUME (UL)	20
REAGENT VOLUME (UL)	100
INCUBATION TIME (SEC)	0
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	30
TIME INTERVAL (SEC)	20
NUMBER OF READINGS	16
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.
Instrument Application

Analyzer: COBAS BIO
Test: Glucose Hexokinase
Catalog # : G7518/G7517

UNITS	MG/DL
CALCULATION FACTOR	0
STANDARD 1 CONC.	*
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	600
TEMPERATURE (DEG C)	37
TYPE OF ANALYSIS	5
WAVELENGTH (NM)	340
SAMPLE VOLUME (UL)	3
DILUENT VOLUME (UL)	30
REAGENT VOLUME (UL)	250
INCUBATION TIME (SEC)	0
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	0.5
TIME INTERVAL (SEC)	150
NUMBER OF READINGS	2
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

Add 13.5 ml distilled water to the 15 ml size vial and 45 ml to the 50 ml size vial.
It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO
Test: Glucose Oxidase
Catalog # : G7519/G7521

UNITS	MG/DL
CALCULATION FACTOR	0
STANDARD 1 CONC.	*
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	500
TEMPERATURE (DEG C)	37
TYPE OF ANALYSIS	5
WAVELENGTH (NM)	500
SAMPLE VOLUME (UL)	4
DILUENT VOLUME (UL)	30
REAGENT VOLUME (UL)	250
INCUBATION TIME (SEC)	0
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	0.5
TIME INTERVAL (SEC)	600
NUMBER OF READINGS	2
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

Add 90 ml distilled water to the 100 ml

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO

Test: Lactate

Catalog # : L7596

UNITS	MG/DL
CALCULATION FACTOR	0
STANDARD 1 CONC.	*
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	User Defined
TEMPERATURE (DEG C)	37
TYPE OF ANALYSIS	6
WAVELENGTH (NM)	546
SAMPLE VOLUME (UL)	6
DILUENT VOLUME (UL)	2
REAGENT VOLUME (UL)	110
INCUBATION TIME (SEC)	30
START REAGENT VOLUME	75
TIME OF FIRST READING (SEC)	0.5
TIME INTERVAL (SEC)	300
NUMBER OF READINGS	1
BLANKING MODE	1
PRINTOUT MODE	2
ALPHA CODE	11

Prepare reagent according to package insert.

Preliminary application

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO
Test: LDH
Catalog # : L7535/I7572

UNITS	U/L
CALCULATION FACTOR	8039
STANDARD 1 CONC.	0
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	0.15
TEMPERATURE (DEG C)	37
TYPE OF ANALYSIS	2
WAVELENGTH (NM)	340
SAMPLE VOLUME (UL)	5
DILUENT VOLUME (UL)	30
REAGENT VOLUME (UL)	200
INCUBATION TIME (SEC)	0
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	50
TIME INTERVAL (SEC)	10
NUMBER OF READINGS	16
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

Add 13 ml distilled water to the 15 ml size vial and 43 ml to the 50 ml size vial.
It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO
Test: Magnesium
Catalog # : M7527

UNITS	mEq/L
CALCULATION FACTOR	0
STANDARD 1 CONC.	*
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	4
TEMPERATURE (DEG C)	30
TYPE OF ANALYSIS	1
WAVELENGTH (NM)	530
SAMPLE VOLUME (UL)	3
DILUENT VOLUME (UL)	10
REAGENT VOLUME (UL)	300
INCUBATION TIME (SEC)	10
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	0.5
TIME INTERVAL (SEC)	180
NUMBER OF READINGS	1
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

Prepare reagent according to package insert.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO

Test: Inorganic Phosphorus

Catalog # : P7516

UNITS	MG/DL
CALCULATION FACTOR	0
STANDARD 1 CONC.	*
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	12
TEMPERATURE (DEG C)	25
TYPE OF ANALYSIS	5
WAVELENGTH (NM)	340
SAMPLE VOLUME (UL)	5
DILUENT VOLUME (UL)	50
REAGENT VOLUME (UL)	250
INCUBATION TIME (SEC)	0
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	0.5
TIME INTERVAL (SEC)	160
NUMBER OF READINGS	2
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

Prepare reagent according to package insert.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO
Test: Total Protein
Catalog # : T7528

UNITS	G/DL
CALCULATION FACTOR	0
STANDARD 1 CONC.	*
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	8.5
TEMPERATURE (DEG C)	25
TYPE OF ANALYSIS	1
WAVELENGTH (NM)	540
SAMPLE VOLUME (UL)	5
DILUENT VOLUME (UL)	40
REAGENT VOLUME (UL)	250
INCUBATION TIME (SEC)	10
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	0.5
TIME INTERVAL (SEC)	180
NUMBER OF READINGS	2
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

Reagent supplied as a ready to use liquid

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO
Test: Triglyceride (GPO)
Catalog # : T7531

UNITS	MG/DL
CALCULATION FACTOR	0
STANDARD 1 CONC.	*
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	1000
TEMPERATURE (DEG C)	37
TYPE OF ANALYSIS	5
WAVELENGTH (NM)	485
SAMPLE VOLUME (UL)	2
DILUENT VOLUME (UL)	30
REAGENT VOLUME (UL)	250
INCUBATION TIME (SEC)	0
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	0.5
TIME INTERVAL (SEC)	300
NUMBER OF READINGS	2
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

Add 13 ml distilled water to the 15 ml size vial and 43 ml to the 50 ml size vial.
It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO
Test: Uric Acid
Catalog # : U7580/U7581

UNITS	MG/DL
CALCULATION FACTOR	0
STANDARD 1 CONC.	*
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	25
TEMPERATURE (DEG C)	37
TYPE OF ANALYSIS	5
WAVELENGTH (NM)	520
SAMPLE VOLUME (UL)	6
DILUENT VOLUME (UL)	30
REAGENT VOLUME (UL)	250
INCUBATION TIME (SEC)	0
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	0.5
TIME INTERVAL (SEC)	300
NUMBER OF READINGS	2
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

Add 13.5 ml distilled water to the 15 ml size vial and 45 ml to the 50 ml size vial.
It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO

Test: Triglyceride

Catalog # : T7532

UNITS	MG/DL
CALCULATION FACTOR	0
STANDARD 1 CONC.	*
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	700
TEMPERATURE (DEG C)	37
TYPE OF ANALYSIS	1
WAVELENGTH (NM)	500
SAMPLE VOLUME (UL)	3
DILUENT VOLUME (UL)	30
REAGENT VOLUME (UL)	250
INCUBATION TIME (SEC)	10
START REAGENT VOLUME	0
TIME OF FIRST READING (SEC)	0.5
TIME INTERVAL (SEC)	600
NUMBER OF READINGS	2
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	11

Reconstitute reagent to 1.5 X the volume stated on the vial label.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.

Instrument Application

Analyzer: COBAS BIO
Test: Iron
Catalog # : I7505/I7506

UNITS	UG/DL
CALCULATION FACTOR	0
STANDARD 1 CONC.	500
STANDARD 2 CONC.	0
STANDARD 3 CONC.	0
LIMIT	500
TEMPERATURE (DEG C)	37
TYPE OF ANALYSIS	6
WAVELENGTH (NM)	560
SAMPLE VOLUME (UL)	50
DILUENT VOLUME (UL)	40
REAGENT VOLUME (UL)	240
INCUBATION TIME (SEC)	60
START REAGENT VOLUME	5
TIME OF FIRST READING (SEC)	5
TIME INTERVAL (SEC)	540
NUMBER OF READINGS	2
BLANKING MODE	1
PRINTOUT MODE	1
ALPHA CODE	Positive

Prepare reagent according to package insert.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

* Insert standard value if required

Rev. 12-02

Pointe Scientific, Inc.
Instrument Application

Analyzer: COBAS BIO
Test: Beta Hydroxybutyrate
Catalog # : H7587

UNITS	MMOL/L
CALCULATION FACTOR	0
STANDARD 1 CONC.	1
STANDARD 2 CONC.	1
STANDARD 3 CONC.	1
LIMIT	0
TEMPERATURE (DEG C)	37
TYPE OF ANALYSIS	6
WAVELENGTH (NM)	505
SAMPLE VOLUME (UL)	6
DILUENT VOLUME (UL)	14
REAGENT VOLUME (UL)	215
INCUBATION TIME (SEC)	300
START REAGENT VOLUME	36
TIME OF FIRST READING (SEC)	0.5
TIME INTERVAL (SEC)	60
NUMBER OF READINGS	5
BLANKING MODE	1
PRINTOUT MODE	3
ALPHA CODE	

It is recommended that two levels of control material be assayed daily.

* Insert standard value if required